

FOOD / CATERING
Detectable Blue Strapping




Ref	Description	Quantity
15682	Detectable Blue Tape 2.5cm x 5m	Spool

Finger Stalls



Ref	Description	Quantity
1150	Finger Stalls Blue ML (6)	Pk

Finger Bobs



Ref	Description	Quantity
13061	Finger Bobs (5) Medium	Pk of 5

FOOD / CATERING
Finger Extension Plasters



Ref	Description	Quantity
1603S	Finger Extension Plasters Fabric	50
17914S	Finger Extension Plasters Wash proof	50
15914S	Finger Extension Plasters Blue	50

Sterochef Blue Detectable Dressings



Ref	Description	Quantity
15640S	Sterochef Dect. Dressings (Asst.) Sizes (100)	Box
15603S	Sterochef Dect. Dressings 7.5cm x 2.5cm (100)	Box
15605S	Sterochef Dect. Dressings 7.5cm x 5cm (50)	Box
15913S	Sterochef Dect. Dressings (Anchor) (50)	Box
15912S	Sterochef Dect. Dressings (F/Tip) (50)	Box
15615/16	Sterochef Dect. Dressings Asst (20)	Box
15914S	Finger Extension Plaster (Blue) (50)	Box

Hypoallergenic adhesive. Complies with food & hygiene regulations (SI1990 No 1431). Metal detectable. Comprehensive range. Sterile. Individually wrapped. Comfortable & washproof.

Burn Free Bottles




Ref	Description	Quantity
BF1	Burn Free Bottle 118 ml	Single
BF2	Burn Free Bottle 48 ml	Single

BURNFREE Pain Relieving Gel is available in 48ml and 118 ml Bottles.

This gel is ideal for all minor burns. It is the same solution found in the burn dressings and fire/trauma blankets. When using BURNFREE Pain Relieving Gel, gently apply a liberal amount to the injured area without rubbing it into the wound. In order for the gel to be fully effective, it is necessary to apply a thick layer over the damaged area.

The dressings are more appropriate for larger and more serious burns. Since the open-cell foam dressings will maintain a thick covering of gel on the burn, they are a more efficient means of covering the area.

FOOD / CATERING
Burn Free Sachets



Ref	Description	Quantity
BF3	Burn Free sachet	Single

Burn Free Dressing



Ref	Description	Quantity
BF4	Burn Free Dressing 10cm x 10cm	Single
BF6	Burn Free Dressing 5cm x 15cm	Single
BF7	Burn Free Dressing 20cm x 20cm	Single
BF8	Burn Free Dressing 30cm x 30cm	Single
BF9	Burn Free Dressing 5cm x 5cm	Single

BURNFREE Sterile Burn Dressings are designed for first aid use on all types of burn injuries. The dressings are constructed of a special, soft, open-cell foam saturated with the BURNFREE solution. BURNFREE dressings have several beneficial physical characteristics such as:

Efficiently covers the burned area and yet elastic enough to accommodate possible swelling or the need to flex joints.
They are light, easily transportable sterile dressings.

A single application of this dressing is an effective first-aid measure on burns.

The dressings have a foam carrier and the foam's high absorption capabilities allow complete delivery of the gel to the burn site. BURNFREE dressings are non-adherent and the inert foam base will not shed threads, strands, fibers, etc., that might become embedded in the wound.

The moist BURNFREE solution facilitates easy removal for inspection and treatment of the injured area. The gel leaves no residue. It can easily be removed with water during normal debridement. BURNFREE gel is completely water soluble.

Burn Free Kit



Ref	Description	Quantity
BK1	Burn Kit - Contents Include: 1 x Burn Kit Box 1 x Paramedic Shears 1 x 120ml Bottle Of Burnfree Cooljel 118ml 2 x 10cm x 10cm Burnfree Dressings 2 x Pks (3) Burnfree Sachets 2 x 3" Conform Bandages 1 x 1" Zinc Oxide Tape 2 x Steropad Non-Adherent Dressings 1 x Pair Of Latex Gloves	Single

FOOD / CATERING
Blue Gloves



Ref	Description	Quantity
	This page is currently being updated	